

Recommended for
educational purposes

MX 05

Economical monocular
microscope

- Triple objective nosepiece
- 45° inclined head
- 3 objectives achromat: 4x/0.10, 10x/0.25, 40x/0.65
- Widefield eyepieces: 10x/16 mm and 25x
- Built-in LED illumination adjustable 5 V, 1 W
- Up and down illumination



Specifications

MX 05 monocular microscope

Standard set

Viewing Head	monocular head, Inclined at 45°
Eyepiece	WF10x/16 and WF25x
Objective	Achromat: 4x/0.10, 10x/0.25, 40x/0.60
Condenser	NA0.65 with Disc Diaphragm
Nosepiece	triple nosepiece
Stage	Plain Stage with slide clips 95x95 mm
Illumination	up and down LED illumination
Focusing System	coaxial coarse and fine adjustment
Additionally (on request)	objectives, digital cameras, software for management of digital albums

Ordering Information

Description

Code

MX 05 monocular microscope, standard set

09.0005.01



We reserve the right to change specification without notice.



West Medica Produktions- und Handels-GmbH
Brown-Boveri-Straße 6, B17-1,
2351 Wiener Neudorf, Austria
tel.: +43 (0) 2236 892465, fax: +43 (0) 2236 892464
vienna@westmedica.com, www.microoptix.com

Official Distributor